

CleanClever Features Diagram

PFA Vessel

- Various sizes for sub-boiling vapor cleaning

PTFE CleanClever Sleeve

- Threaded union connects vessel to bottom reservoir
- Customized sleeves for various sized vessels & bottom reservoirs

PFA Bottom Reservoir

- Holds acid or water used for vapor cleaning

CleanClever Well

- Heated well that holds PFA bottom reservoir

Heated Surface

- FluoroPlate
- FluoroPlate with EasyDigest Rack
- Wells on top of EvapoClean & CleanAcids

