

DX Series Autosamplers

DXCi Autosamplers

Parabolic Probe

- Robust, heavy-walled, self-supporting

Position Autocorrection

- Automatically goes to the correct position every time
- Automatically retries position if an obstacle is encountered
- Automatically reports unmoved obstruction

Highest Speed of Operation

- Corrects for stall events, even after years of use
- Keeps analysis running, even after obstruction

DXCi / DX Autosamplers

FAST Ready

- FAST sample loading system installed

Dual Independent Rinse Pumps

- Each station can use a different rinse solution

DXCi / DX / DXe Autosamplers

Upgrade Options

- SampleSense FAST for ultra-high sample throughput
- prepFAST X for automatic ICP / ICPMS calibrations and automatic inline sample dilution
- Fluoronetic drive for highest chemical resistance

AutoAlign Arm

- Self-realigns after encountering capped sample tube
- Metal-free fluoropolymer probe aligner

Sample Containers

- Flexible rack configurations
- Samples from micro vials to large bottles

Dual Overflowing Rinse Stations

- Ultra-low cross contamination
- Optional organic solvent rinse
- Low rinse consumption

Chemically Resistant

- Designed for exceptionally long lifetime
- No metallic components over samples

